

## STANDARD OPERATING PROCEDURE

**Title: Immunofluorescence (IF) Microscopy using EVOS™ M5000 Imaging System**

**SOP#: M-137**

**Version #: 2**

**Author: R. Roberts**

**Date Approved: September 24, 2020**

**Date Modified: January 1, 2023**

### 1. PURPOSE

This procedure will describe how to prepare antibodies/protein targets for immunofluorescence (IF). The procedure utilizes fluorescent-labeled antibodies to detect specific target proteins at the subcellular level, followed by imaging using the EVOS™ M5000 Imaging System to provide information such as protein expression and localization.

### 2. SCOPE

This procedure applies to all samples that are processed for immunofluorescence analysis using EVOS™ M5000 Imaging System. This procedure incorporates the Human Protein Atlas (HPA) IF staining protocol (with minor changes as described in this SOP) for adherent and suspension cells using a 96-well microplate format but uses the EVOS™ M5000 Imaging System. The fixation and permeabilization of cells, primary and secondary antibody incubations 4',6-diamidino-2-phenylindole (DAPI) fluorescent nuclear staining and mounting procedures will be described in this SOP.

### 3. RESPONSIBILITIES

It is the responsibility of person(s) performing this procedure to be familiar with cell culture and paraformaldehyde (PFA) lab safety procedures and the use of the EVOS™ M5000 imaging system. The interpretation of results must be done by a person trained in the procedure and familiar with such interpretation. It is the responsibility of the analyst to follow the procedure and document any deviations and all observations in a laboratory notebook.

#### 4. EQUIPMENT

- Invitrogen™ EVOS™ M5000 Imaging System; Thermo Fisher Scientific
- Standard Laboratory Cell Culture CO<sub>2</sub> Incubator
- Biosafety Cabinet/Hood Class II
- Laboratory Chemical Fume Hood
- Centrifuge with swing rotor for tubes and microplates, Eppendorf, Model #5810 R
- Vortex
- pH meter
- Stirring Hot plate

#### 5. MATERIALS

Use the following recommended materials, or equivalent:

- Greiner Sensoplate™ glass bottom multi-well plates; Millipore Sigma, Cat. #M4187-16EA
- Shi-fix™ coated 96-well microplate, Everest Biotech, Cat. #SB-Shifix96
- 250 ml beaker and magnetic stir bars
- Concentrated NaOH
- Concentrated HCl
- Deionized (DI) Water (MilliQ or better)
- 10X Phosphate Buffered Saline (PBS), Fisher Scientific, Cat. #BP399-1 (diluted to 1X with DI water)
- Superfibronectin; Millipore Sigma, Cat. #S5171-.5MG
- Heat Inactivated Fetal Bovine Serum (FBS), Life Technologies, Cat. #16140071.
- 16% Paraformaldehyde (PFA); VWR, Cat. #BT140770-10X10
  - **Prepare 4% PFA/1X PBS with FBS as follows:**  
*Always use lab coat, eye protection and gloves and perform work in a fume hood when preparing and using 4% PFA /1X PBS with FBS.*
    1. Pre-heat 75 ml 10% FBS diluted with 1X PBS in a 250 ml beaker to 60°C on a hot plate with magnetic stirrer.
    2. Add 25 ml 16% PFA to the warm FBS/PBS (Step 1) and continue stirring for 20 min at 60°C.
    3. Slowly add concentrated NaOH to the solution from Step 2 until the solution reaches pH ~11.

4. Allow the solution to cool down to room temperature (~20°C).
  5. Adjust pH to 7.2 - 7.3 using first concentrated HCl, and then diluted HCl for fine adjustment.
  6. Aliquot solution and store at -20°C until use. Thaw at room temperature right before use.
- Triton X-100 Solution; Millipore Sigma, Cat. # 93443-100ML.
    - **Prepare at time of use - Permeabilization buffer solution as follows:**
      - Mix 0.05 mL of TritonX-100 solution with 4.95 mL of 1X PBS
    - **Prepare at time of use - Blocking buffer solution as follows:**
      - Mix 0.25 mL of FBS with 2.25 mL of Permeabilization buffer
  - DAPI; VWR, Cat. #10180-614.
    - **Prepare DAPI solution as follows:**
      - Dissolve 10 mg of DAPI powder in 2 ml DI water to a 5 mg/ml (14.3 mM) concentrated stock solution. Prepare 100 uL aliquots and store at -20°C until needed.
      - When ready to use thaw 5 mg/ml DAPI stock solution at room temperature, then combine 80 µl of 5 mg/ml DAPI stock solution with 920 µl DI water to a pre-diluted 400 µg/ml DAPI solution. Discard any unused 5 mg/ml DAPI stock solution. Prepare 5 uL aliquots of the DAPI solution (400 ug/ml) and store at -20°C until needed.
      - When ready to use thaw pre-diluted 400 µg/ml DAPI solution at room temperature to prepare at time of use, 0.2 µg/ml DAPI with 1X PBS (0.2 µg/ml DAPI/1X PBS) for experiment. The 400 µg/ml DAPI solution can be refrozen up to 2 times, then it should be discarded.
  - Aluminum PCR Plate Sealers, Greiner Bio-One; VWR, Cat. # 82050-998
  - Pipettes and tips
  - Micro centrifuge and tubes
  - 50 ml 0.2 µm filter units; Thermo Fisher Scientific, Cat. #564-0020
  - Glycerol; Thermo Fisher Scientific, Cat. #15514011
    - **Prepare Glycerol/10X PBS as follows:**
      - Add 5 ml 10X PBS to 45 ml of glycerol. Mix by inverting.
      - Autoclave or filter sterilize.
      - Store at Room Temperature

## 6. **REAGENTS**

### 6.1. **For Rabbit or Mouse Monoclonal Antibody Evaluation**

- Cells corresponding to the antibody or protein of interest such as cell lines (suspension or adherent) or equivalent.
- Purified monoclonal antibody corresponding to the protein of interest
- Mouse or rabbit anti-alpha tubulin marker; Abcam, Cat. #Ab7291 or Ab176560, respectively.
- Goat anti-rabbit Alexa Fluor 488 and goat anti-mouse Alexa Fluor 555 secondary antibodies; Thermo Fisher Scientific, Cat. #A11034 and A21424, respectively or other species as appropriate conjugated Alexa Fluor secondary antibodies.

## 7. **PROCEDURE**

### 7.1. **Guidelines**

- All volumes in the protocol refer to volume per well in a 96-well microplate. Unspecified incubation temperatures are at room temperature. Immunostaining can be performed manually or using a liquid handling robot.

### 7.2. **Preparation of Primary Antibodies**

- Before immunostaining, dilute primary antibodies to 1:150 in Blocking Buffer (10% FBS constituted in 0.1% Triton X-100/1X PBS) together with 1:1000 dilution anti-alpha tubulin marker in one tube.

### 7.3. **Preparation of Secondary Antibodies**

- Before immunostaining, dilute secondary antibodies to 1:300 in Blocking Buffer (10% FBS constituted in 0.1% Triton X-100/1X PBS) together with 0.2 ug/ml DAPI/1X PBS in one tube and protect from light.

### 7.4. **Cell Culture**

- Grow adherent cells overnight in 100 µl cell culture media (according to cell provider's instructions) on fibronectin (12.5 µg/ml) coated 96-well glass bottom microplate (precoated according to manufacturer instructions). Note, optimal seeding concentration/confluency for IF staining varies between cell lines and test samples it may be required to define optimal conditions for desired cell line/target. For suspension cells (10 million cells/mL in

1X PBS) use the Shi-fix™ coated 96-well microplate at a recommended concentration between 200,000 -500,000 cells/well. Follow Shi-fix™ manufacturer's instructions and proceed to fixation. Note, optimal cell count for IF staining varies between cell lines and test samples it may be required to define optimal conditions for desired cell line/target.

#### **7.5. Fixation**

- Gently remove cell culture media (aspirate ~ 100 µl from all wells).
- Fix cells by incubating with 40 µl 4% PFA/10% FBS/1X PBS for 15 minutes.
- Gently wash cells with 40 µl 1X PBS for 5 minutes repeat 3X.

#### **7.6. Permeabilization**

- Add the Permeabilize solution to the cells (40 µl 0.1% Triton X-100/1X PBS) and incubate for 15 minutes.

#### **7.7. Blocking**

- Gently remove Permeabilization solution from wells.
- Add 40 µl Blocking Buffer (10% FBS constituted in 0.1% Triton X-100/1X PBS) to the cells and incubate for 1 hour.

#### **7.8. Primary Antibody incubation**

- Gently remove Blocking Buffer from wells.
- Add 40 µl Blocking Buffer (10% FBS constituted in 0.1% Triton X-100/1X PBS) with diluted primary antibodies and anti-alpha tubulin marker (Prepared in Step 7.2).
- Incubate overnight at 4°C.

#### **7.9. Secondary Antibody and DAPI Nuclear Staining incubation**

- Gently remove Blocking Buffer containing primary antibodies and gently wash cells with 40 µl 1X PBS for 5 minutes repeat 3X.
- Add 40 µl Blocking Buffer containing diluted secondary antibodies and DAPI (Prepared in Step 7.3).
- Incubate for 1 hour in the dark.

#### **7.10. Mounting**

- Gently remove Blocking Buffer containing secondary antibodies and DAPI.
- Gently wash cells with 40 µl 1X PBS for 5 minutes repeat 3X.

- Gently add glycerol/10X PBS to cover the cells (~240  $\mu$ l) and seal the plate with an adhesive aluminum PCR plate seal to store if not immediately imaging.
- Image the cells immediately or store plates at 4°C for no more than two weeks before imaging for best results.

## 8. **REFERENCED DOCUMENTS**

- 8.1 HPA Standard Immunostaining Protocol, version Feb 2019  
<https://www.protocols.io/view/hpa-cell-atlas-standard-immunostaining-protocol-x2dfqa6>
- 8.2 Shi-fix™ coated 96-well microplate protocol, Everest Biotech Product Datasheet