

Antibody Pairing Results for CPTC-ANXA1

“Capture” (bottom) Antibody	“Indicator” (top) Antibody		
	CPTC-ANXA1-1	CPTC-ANXA1-2	CPTC-ANXA1-3
CPTC-ANXA1-1	-	-	+++
CPTC-ANXA1-2	+++	-	+
CPTC-ANXA1-3	+++	-	-

Method – BioRad ProteOn XPR:

- Surface coated with Rabbit anti-mouse IgG Fc specific antibody
- “Bottom” antibody flowed over the surface
- Antigen flowed over the surface
- Orthogonal flow of “Top” antibody over the surface

Results:

- No pairing observed
- + Weak Pairing
- ++ Moderate Pairing
- +++ Strong pairing